

## Best Localization Needles



### Best Medical International supplies breast localization needles with Iodine-125 seeds

- Best localization needles are available in 5, 7, 10 and 15 cm lengths
- Available in various ranges of seed activities with a shelf life of 90 days
- Needles are supplied sterile and are individually packaged and shielded

© 2015 Best Medical International, Inc.


Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA  
703 451 2378 800 336 4970 [www.bestmedical.com](http://www.bestmedical.com) [www.teambest.com](http://www.teambest.com)

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA




## Best Double-Balloon Breast Brachytherapy Applicator

The Best Double-Balloon Breast Brachytherapy Applicator utilizes an outer balloon to expand within and shape the resection cavity, and an inner balloon which expands independently of the outer balloon. This allows physicians to position 4 treatment catheters within the outer balloon, optimizing both dose conformity and homogeneity.



- Improved dose distribution and conformity
- Less dose to critical organs such as skin, lung, heart, chest wall, etc.
- Convenient to use
- US FDA 510(k) registered



© 2015 Best Medical International, Inc.

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA  
703 451 2378 800 336 4970 [www.bestmedical.com](http://www.bestmedical.com) [www.teambest.com](http://www.teambest.com)

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA

