

Best™ *medical international*

Best™ Total Solutions for Breast Cancer Treatment

TEAMBEST GLOBAL



www.teambest.com

Best™ HDR Afterloader

- Provides 20 channels for dose delivery
- Battery back-up in case of power failure and automated wire recovery
- Source wire tracking using sensors and encoders
- Quick source replacement process — reducing down-time
- Available with Cobalt-60 or Iridium-192 sources. Also, Ytterbium-169 source available soon!





Best™ E-Beam™ Robotic IORT Linac System

- Mobile intraoperative radiation therapy for use in operating rooms
- Targets residual cancer cells during surgery, reducing the need for radiotherapy or chemotherapy at a later date in most cases

Best™ Breast Double-Balloon Brachytherapy Applicator

- Improved dose distribution and conformality
- Less dose to critical organs such as skin, lung, heart, chest wall, etc.
- Easy to use and available in 3 different sizes
- US FDA 510(k) registered



Best Medical International, Inc. © 2020

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
tel: 703 451 2378 800 336 4970 www.besttotalsolutions.com www.teambest.com

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA

* Certain products shown are not available for sale currently.