


Best™ ABT Molecular Imaging

BG-75 BIOMARKER GENERATOR



The BG-75 Biomarker Generator integrates a compact mini-cyclotron, kit based micro-chemistry, and automated quality control, simplifying in-house production of ^{18}F -FDG and advanced biomarkers.

- Push button graphic interface
- Kit based chemistry
- Single or batch dose production
- Final dose delivery to syringe or vial (option)
- Automated quality control testing
- Integrated cyclotron & chemistry self-shielding
- Complete production lab in a 30 x 30 meter area



Best ABT Molecular Imaging, Inc. © 2019

Best ABT Molecular Imaging 3024 Topside Business Park Drive, Louisville, TN 37777
tel: 865 982 0098 800 336 4970 www.bestabt.com www.teambest.com

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA

