



BestTM ABT *Molecular Imaging*



- **BG-75 Biomarker Generator is a compact cyclotron**
- **Shelf-shielded, integrated unit for a small room (5 x 5 meters)**
- **Easy-to-operate Push Button System to deliver “Dose-On-Demand”**
- **On-board chemistry and QC units for in house production of F(18)-FDG and other PET agents**



tel: 865.982.0098 fax: 865.982.0411
www.bestabt.com www.teambest.com

TEAMBEST GLOBAL



www.teambest.com