


NEW! Best Model B35adp Alpha/Deuteron/Proton Cyclotron for Medical Radioisotope Production and Other Applications

- Proton Particle Beam: 30 MeV Energy, 500 μ A Beam Current
- Deuteron Particle Beam: 15 MeV Energy, 250 μ A Beam Current
- Alpha Particle Beam: 35 MeV Energy, 100 μ A Beam Current




NEW! Best Model 150p Cyclotron for Proton Therapy (Patent Pending)

- From 70 MeV up to 150 MeV Non-Variable Energy
- Dedicated for Proton Therapy with two beam lines & two treatment rooms
- For all Medical Treatments including: Benign & Malignant Tumors for Neurological, Eye, Head/Neck, Pediatric, Lung Cancers, Vascular/Cardiac/Stenosis/Ablation, etc.



Best Particle Therapy 400 MeV ion Rapid Cycling Medical Synchrotron (iRCMS) for Proton-to-Carbon, Variable Energy Heavy Ion Therapy, with or without Gantries — Upgradeable Single and Multi-Room Solutions



Upgradeable Single Room & Multi-Room Solutions

- Intrinsically small beams facilitating beam delivery with precision
- Small beam sizes – small magnets, light gantries – smaller footprint
- Highly efficient single turn extraction
- Efficient extraction – less shielding
- Flexibility – heavy ion beam therapy (protons and/or carbon), beam delivery modalities



BEST RADIATION THERAPY & DIAGNOSTIC CENTER

ion Rapid Cycling Medical Synchrotron (iRCMS) Racetrack

www.bestcyclotron.com • www.bestabt.com • www.bestproton.com
 BCS tel: 604 681 3327 • Best ABT tel: 865 982 0098 • BPT tel: 703 451 2378

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA