







Best[®] medical international

Applicators



healthcare for everyone
TeamBest[®]
 Your True Partner

**Best Medical can
provide any standard
or customized
configuration made
to your exact
specifications, sterile
or non-sterile,
shipped within 24
hours — worldwide!**

Contact us:

**Best Medical International
7643 Fullerton Road
Springfield, Virginia 22153 USA**

Tel: 703-451-2378 • 800-336-4970

Fax: 703-451-5228

www.bestmedical.com • www.teambest.com



Mission Statement

Everyone deserves the best health care. Our goal is to work with medical professionals to provide the Best products and services. As a leader in innovative technology, our mission is to uphold our reputation for excellence in the health care field by developing, manufacturing, and delivering cost-effective, high quality products to benefit patients around the world.

Corporate Values

Fulfilling our mission requires a dedicated employee team—one that is committed to uncompromised quality and safety, and fostering a sense of pride in the company.

Best employees will:

- *Honor the dignity and uniqueness of each person, including fellow employees, customers, distributors, and vendors, through their actions, words, and behavior. They will treat others with respect, as they would wish to be treated.*
- *Embrace the spirit of cooperation and teamwork for which Best Medical is known, an environment which supports and recognizes effort and results.*
- *Promote and practice open and honest communication at all times.*
- *Listen to others and remain open to suggestions at all times, as we focus on achieving excellence in our work.*
- *Attract and retain leaders with exceptional abilities and integrity.*
- *Ensure all resources are managed in a responsible and effective manner.*

Best ADVANTAGE

- Best Medical International is the **ONLY** supplier of Iridium-192 and Gold-198 seeds.
- We are the first and only provider of Iridium-192 in nylon ribbons in a wide range of activities.
- We are the first and only company to provide custom service for Gold-198 permanent implants, as well as a wide range of activities.
- We are the first and only company that will call our customers the morning after shipping to inform them of mode of transportation and expected arrival time.
- We can supply Ir-192, Au-198, Pd-103, I-125 and Brachytherapy accessories on a single Purchase Order.
- We ship most materials within 24 hours of receipt of the order.

Catalog 2013 Table of Contents

- Round-Handle Fletcher-Suit-Delclos Applicator Set (640R) 1
- Round Handle Fletcher-Suit-Delclos Applicator Set with Bracket (640RB) 2
- Fletcher-Suit-Delclos Mini-Ovoid Applicator Set (640M) 3
- Delclos Mini-Ovoid Applicator Set with Bracket Assembly (640MRB) 4
- Burnett Afterload Vaginal Applicator (700A) 5
- Delclos Uterine/Vaginal Afterloading Applicator and Dome Assembly (670) 6
- Dummy Sources (Cesium) 7
- Mick® 200-TPV Applicator Kit 8
- Ordering Information: Optional Instrumentation and Accessories 9




Round-Handle Fletcher-Suit-Delclos Applicator Set (640R)

Part # 5005-0082

Set Contains:

3	(a) Central Tandems
3	(b) Cervical Stops (1 with flag, 2 without flag)
3	(c) Plastic Tandems and Pushers
1	(d) Storage/Sterilization Tray (not shown)
1	(e) Pair Colpostat w
2	(f) Bucket Type Source Carriers
1	(g) Pair Nylon Caps 2.5 cm Diameter for Ovoids
1	(h) Pair Nylon Caps 3.0 cm Diameter for Ovoids
1	(i) Allen Wrench




CAUTION: This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.

Round Handle Fletcher-Suit-Delclos Applicator Set with Bracket (640RB)

Part # 5005-0083

Set Contains:

3	(a) Central Tandems
3	(b) Cervical Stops (1 with flag, 2 without flag)
3	(c) Plastic Tandems and Pushers
2	(d) Allen Wrenches
1	(e) Storage/Sterilization Tray (not shown)
1	(f) Pair Colpostat with Bracket (to hold central tandem along with colpostat and also to separate the ovoids)
1	(g) Pair Bucket Type Source Carriers
1	(h) Pair Nylon Caps 2.5 cm Diameter for Ovoids
1	(i) Pair Nylon Caps 3.0 cm Diameter for Ovoids




CAUTION: This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.

Fletcher-Suit-Delclos Mini-Ovoid Applicator Set (640M)

Part # 5005-0084

Set Contains:

3	(a) Central Tandems
3	(b) Cervical Stops (1 with flag, 2 without flag)
3	(c) Plastic Tandems and Pushers
1	(d) Allen Wrench
1	(e) Storage/Sterilization Tray (not shown)
1	(f) Pair Colpostat
2	(g) Bucket Type Source Carriers



* Optional: Shielded Caps 2 cm, 2.5 cm, 3 cm Diameters.


CAUTION: This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.

Delclos Mini-Ovoid Applicator Set with Bracket Assembly (640MRB)

Part # 5005-0085

Set Contains:

3	(a) Central Tandems
3	(b) Cervical Stops (1 with flag, 2 without flag)
3	(c) Plastic Tandems and Pushers
2	(d) Allen Wrenches
1	(e) Storage/Sterilization Tray (not shown)
1	(f) Pair Colpostat with Bracket
2	(g) Bucket Type Source Carriers



*** Optional: Shielded Caps 2 cm, 2.5 cm, 3 cm Diameters.**


CAUTION: This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.

Burnett Afterload Vaginal Applicator (700A)

Part # 5005-0088

Set Contains:

8	Cylon Cylinders
4	Plastic tandems with metal screws and special long handle (cap holder) and tube source adapter assembly arranged in a stainless steel tray.



External Dimensions of Cylinders:

A	20 x 66 mm
B	20 x 87 mm
C	25 x 87 mm
D	30 x 87 mm
E	35 x 87 mm
F	25 x 95 mm
G	30 x 95 mm
H	35 x 95 mm


CAUTION: This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.

Delclos Uterine/Vaginal Afterloading Applicator and Dome Assembly (670)

Part # 5005-0089

Set Contains:

6	Sets of Cylinders with Diameters 20 mm, 25 mm, 30 mm, 35 mm, 40 mm and 45 mm
6	Domes (matching cylinder diameters)
3	Uterine Tandems
1	Vaginal Tandem (dome assembly)
1	Vaginal Tandem (with a flat end)
1	Storage/Sterilization Tray
1	Standard Allen Wrench
1	Extended Allen Wrench
1	Cylinder Extractor
1	Cervical Stop
4	Plastic Tandems and Pushers




CAUTION: This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.

Dummy Sources (Cesium)

Part # 513-1

These are Dummy Tungsten Sources with a length of 20 mm and a diameter of 3 mm to simulate Cesium sources, arranged with 20 mm long spacers in a plastic tube so as to load conveniently in a central tandem. Because of the high density of Tungsten a better Radiograph is obtained. Custom arrangements of these sources are available.

Four (4) Cesium Dummy Sources Arranged in a Plastic Tube



Dummy Sources (Cesium)

Part # 513-2

Instead of total length of 20 mm, active length of 14 mm also could be used for loading in plastic tube. Custom arrangements of these are available.

Also Available:




513	Cesium Dummy Sources for 3M Sources
696	Cesium Dummy Sources for Amersham Sources

Other Custom Loadings Available



CAUTION: This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.

Mick[®] 200-TPV Applicator Kit		
Part # 5005-9707		
Set Includes:		Part #
1 ea.	Mick 200-TPV Applicator (instrument only)	5005-9701-1
1 ea.	8x Mick TP Needles, 17 or 18 Gauge "Authentic" (carton)	10203-25 (17G) or 10209-25 (18G)
1 ea.	Diddler (needle adjustment tool)	5002-0010
1 ea.	Seed Loading V-Block with Seed Gauge	5005-7509
1 ea.	Pushwire Rod Assembly (spare)	5005-9707-PW
2 ea.	Reverse Action Tweezer	5002-0009
1 ea.	Loop Wire, 22 cm (needle cleaning device)	5005-7702
1 ea.	Stainless Steel Screwdriver/Allen Wrench Combo	5005-7311-C
1 ea.	Sterilization & Storage Cassette	5005-0505
1 ea.	Instruction Manual	405-06



Steerable "Authentic" Mick TP Needles!

Circumvent pubic arch interference with the use of the Mick Needle Former.

Ordering Information: Optional Instrumentation and Accessories

Part # 5005-9707

Product:	Part #
Reverse Action Tweezer	5002-0009
Diddler (needle adjustment tool)	5002-0010
Transperineal Template – Transparent, 17 Gauge	5005-8813-17
Transperineal Template – Transparent, 18 Gauge	5005-8813-18
Mick Needle Former (steerable Mick needles)	5005-0226
V-Block (cartridge seed-loading device)	5005-7509
Stainless Steel Sterilization Tray with Instrument Holder	5005-7308-STI
Wooden Carrying Case	5005-7611
Morgenstern, Stabilizing Needle, 18 Gauge	10207
Practice Phantom with 100 Dummy Seeds	5005-8113-PIP
KGM Passer (radiation shielding device)	5005-9512
Inspector (radiation survey meter)	5005-963ERAD
Springs (for magazines)	5005-9609-SP
E-Rings (for magazines)	5005-9609-E
Prostate Phantom, CIRS	5005-9603-53
Mick Magazine, Reusable (10 seed capacity)	5005-7609
Mick Magazine, Disposable (15 seed capacity)	10113-7
Mick TP Needle, 18 Gauge “Authentic”	10209-25
Mick TP Needle, 17 Gauge “Authentic”	10203-25
Pushwire Rod Assembly (Mick 200-TPV applicator)	5005-9707-PRA
Stainless Steel Allen Wrench/Screwdriver (combination)	5005-7311-C
Brachystand 6 Ultra for Use with E&K or Siemens Stepper	5005-1000
Sure-Point I (floor mount with Siemens stepper)	5005-FMS
Sure-Point II (table mount with Siemens stepper)	5005-UTMS
Disposable Templates, B&K	1830
Leaded Gloves – Radiation Attenuation Gloves	5005-5327005
Ultrasonic Cleaner with Timer and Heater	5005-7908-10TH
Stainless Steel Dummy Seeds	5005-8212-DS

CAUTION: This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.



Best Medical International, Inc. as author/distributor, believes the technical information and statements contained in this catalog to be reliable and correct at the time of publication. However, Best Medical International, Inc., does not assume any liability for statements, information or recommendations based on/or made in reliance on the information contained in this publication.

Data provided herein are for informational purposes only and are not intended for clinical use of any nature. Use of these data in such a manner is the sole responsibility of the user, and the author/distributor of the information makes no claim or inference as to its clinical suitability and/or accuracy. Best Medical International, Inc. assumes no warranty or liability for use of this information or products. Designs and specifications are subject to modification and/or change(s).

© 2013 Best Medical International, Inc.

Terms and Conditions

Ordering: Orders may be placed with a purchase order number and transmitted by FAX, phone or e-mail.


Payment: Within 30 days.

Freight: FOB origin, freight prepaid and charged back through invoice.

Returns: Unused accessories may be returned to Best Medical International, Inc. within 15 days if found unsuitable for a particular application. No return or authorization number is necessary. A note stating the reason for return, name, and address of hospital will expedite processing of any credit.

Note: Best Medical International, Inc. assumes no liability for consequential damages of any kind for the material when used interchangeably with products of other manufacturers/suppliers or for any direct or indirect results and consequences of its use or misuse by the customer.

Best Medical International, Inc., reserves the right to modify the dimensions of accessories/components as well as the prices.



Are you Searching for the Greatest Healthcare Company in the World?



Look to TeamBest®!

World population is increasing by more than 100 million per year. To meet the needs of our aging and growing population, TeamBest® companies are committed to making quality healthcare affordable and accessible globally.

Ready to Serve Globally for All Your Needs

- Arplay Medical
- Best Cyclotron Systems
- Best Medical Africa
- Best Medical Asia
- Best Medical Canada
- Best Medical China
- Best Medical Engineering
- Best Medical Europe
- Best Medical International
- Best Medical Latin America
- Best Medical Middle East
- Best Medical Parts
- Best Medical Research
- Best Medical Russia
- Best NDT
- Best NOMOS
- Best Particle Therapy
- Best Theratronics
- Best Vascular (Novoste)
- CNMC Company
- Genetherapy Pharmaceuticals
- Huestis Diagnostics
- Huestis Medical and more...

For more information, visit www.teambest.com

© 2013 Best Medical International, Inc.

