

# Best<sup>®</sup> medical *international*

## Flexi-Needles



**Best Medical can  
provide any standard  
or customized  
configuration made  
to your exact  
specifications, sterile  
or non-sterile,  
shipped within 24  
hours — worldwide!**

## **Contact us:**

**Best Medical International  
7643 Fullerton Road  
Springfield, Virginia 22153 USA**

**Tel: 703-451-2378 • 800-336-4970**

**Fax: 703-451-5228**

**[www.bestmedical.com](http://www.bestmedical.com) • [www.teambest.com](http://www.teambest.com)**



## Mission Statement

*Everyone deserves the best health care. Our goal is to work with medical professionals to provide the Best products and services. As a leader in innovative technology, our mission is to uphold our reputation for excellence in the health care field by developing, manufacturing, and delivering cost-effective, high quality products to benefit patients around the world.*

## Corporate Values

*Fulfilling our mission requires a dedicated employee team—one that is committed to uncompromised quality and safety, and fostering a sense of pride in the company.*

## Best employees will:

- *Honor the dignity and uniqueness of each person, including fellow employees, customers, distributors, and vendors, through their actions, words, and behavior. They will treat others with respect, as they would wish to be treated.*
- *Embrace the spirit of cooperation and teamwork for which Best Medical is known, an environment which supports and recognizes effort and results.*
- *Promote and practice open and honest communication at all times.*
- *Listen to others and remain open to suggestions at all times, as we focus on achieving excellence in our work.*
- *Attract and retain leaders with exceptional abilities and integrity.*
- *Ensure all resources are managed in a responsible and effective manner.*

# Best ADVANTAGE

- Best Medical International is the ONLY supplier of Iridium-192 and Gold-198 seeds.
- We are the first and only provider of Iridium-192 in nylon ribbons in a wide range of activities.
- We are the first and only company to provide custom service for Gold-198 permanent implants, as well as a wide range of activities.
- We are the first and only company that will call our customers the morning after shipping to inform them of mode of transportation and expected arrival time.
- We can supply Ir-192, Au-198, Pd-103, I-125 and Brachytherapy accessories on a single Purchase Order.
- We ship most materials within 24 hours of receipt of the order.

---

# Catalog 2013 **Table of Contents**

Flexi-Needles ..... 1

Flexi-Needle Combinations ..... 2

Sterile Flexi-Needles ..... 7





## Flexi-Needles



### Description

Flexi-Needles are sharp or blunt, with closed-ended needle guides, and have a nominal 0.3 cm long tip. They are made of biologically compatible material. (Custom lengths are available upon request.)

### Intended Use

Flexi-Needles are inserted transcutaneously, using special obturators. They are designed to be afterloaded with Iridium-192 or Iodine-125 ribbon sources. Two types of caps are available (Luer Cap or Slip Cap) to place over the loaded Flexi-Needle to hold the ribbons in place during irradiation. Simulator sources can be used for localization; obturators are left in place until the time each is afterloaded with the appropriate radioactive source.

### PRECAUTIONS

- CAUTION:** Device is shipped Non-Sterile or Sterile.
- CAUTION:** Flexi-Needles and Obturators are for single use only.
- CAUTION:** Federal law restricts this device to sale by or on the order of a physician.
- CAUTION:** Ensure Obturator is all the way to the tip of the Flexi-Needle prior to insertion.
- CAUTION:** Luer hubs may separate from the Flexi-Needle and Obturators if excessive force is applied.

**WARNING:** Obturators are made of Tungsten or Stainless Steel. Both obturator and Flexi-Needle may break if bone, hard skin, hard surfaces, or excessive forces are encountered.

<b>Flexi-Needle Combinations</b>			
<i>Gauge</i>	<i>Part #</i>	<i>Outer Diameter</i>	<i>Inner Diameter</i>
<b>1. Flexi-Needle Sharp without Luer Lock, without Obturator**</b>			
13 Gauge	251-**	0.095"	0.078"
14 Gauge	252-**	0.083"	0.062"
15 Gauge	253-**	0.078"	0.057"
15.5 Gauge	491-**	0.072"	0.054"
16 Gauge	254-**	0.066"	0.050"
17 Gauge	255-**	0.058"	0.040"
18 Gauge	256-**	0.050"	0.031"
<b>2. Flexi-Needle Sharp, with Luer Lock, without Obturator**</b>			
13 Gauge	263-**	0.095"	0.078"
14 Gauge	264-**	0.083"	0.062"
15 Gauge	265-**	0.078"	0.057"
15.5 Gauge	554-**	0.072"	0.054"
16 Gauge	266-**	0.066"	0.050"
17 Gauge	267-**	0.058"	0.040"
18 Gauge	268-**	0.050"	0.031"
<b>3. Flexi-Needle Sharp, with Luer Lock and Tungsten Obturator**</b>			
14 Gauge	282-**	0.083"	0.062"
15 Gauge	283-**	0.078"	0.057"
15.5 Gauge	555-**	0.072"	0.054"
16 Gauge	284-**	0.066"	0.050"
17 Gauge	285-**	0.058"	0.040"
18 Gauge	286-**	0.050"	0.031"
<b>4. Flexi-Needle Sharp, with Luer Lock and Stainless Steel Obturator**</b>			
13 Gauge	281-**	0.095"	0.078"
14 Gauge	755-**	0.083"	0.062"
15 Gauge	733-**	0.078"	0.057"
15.5 Gauge	1568-**	0.072"	0.054"
16 Gauge	749-**	0.066"	0.050"
17 Gauge	741-**	0.058"	0.040"
18 Gauge	763-**	0.050"	0.031"

**\*\* please specify length**

*continues on next page*

<b>Gauge</b>	<b>Part #</b>	<b>Outer Diameter</b>	<b>Inner Diameter</b>
<b>5. Flexi-Needle Sharp, without Luer Lock, with Tungsten Obturator**</b>			
14 Gauge	355-**	0.083"	0.062"
15 Gauge	356-**	0.078"	0.057"
15.5 Gauge	1136-**	0.072"	0.054"
16 Gauge	357-**	0.066"	0.050"
17 Gauge	358-**	0.058"	0.040"
18 Gauge	359-**	0.050"	0.031"
<b>6. Flexi-Needle Sharp, without Luer Lock, with Stainless Steel Obturator**</b>			
13 Gauge	354-**	0.095"	0.078"
14 Gauge	757-**	0.083"	0.062"
15 Gauge	735-**	0.078"	0.057"
15.5 Gauge	971-**	0.072"	0.054"
16 Gauge	751-**	0.066"	0.050"
17 Gauge	743-**	0.058"	0.040"
18 Gauge	765-**	0.050"	0.031"
<b>7. Flexi-Needle Sharp, without Luer Lock, with Rod (Tungsten Obturator without Luer Lock)**</b>			
14 Gauge	403-**	0.083"	0.062"
15 Gauge	404-**	0.078"	0.057"
16 Gauge	405-**	0.066"	0.050"
17 Gauge	406-**	0.058"	0.040"
18 Gauge	407-**	0.050"	0.031"
<b>8. Flexi-Needle Sharp, without Luer Lock, with Rod (Stainless Steel Obturator without Luer Lock)**</b>			
13 Gauge	402-**	0.095"	0.078"
14 Gauge	759-**	0.083"	0.062"
15 Gauge	737-**	0.078"	0.057"
16 Gauge	753-**	0.066"	0.050"
17 Gauge	745-**	0.058"	0.040"
18 Gauge	767-**	0.050"	0.031"

**\*\* please specify length**

*continues on next page*

Gauge	Part #	Outer Diameter	Inner Diameter
<b>9. Flexi-Needle Blunt End, without Luer Lock, and without Obturator**</b>			
13 Gauge	257-**	0.095"	0.078"
14 Gauge	258-**	0.083"	0.062"
15 Gauge	259-**	0.078"	0.057"
16 Gauge	260-**	0.066"	0.050"
17 Gauge	261-**	0.058"	0.040"
18 Gauge	262-**	0.050"	0.031"
<b>10. Flexi-Needle Blunt End with Luer Lock, without Obturator**</b>			
13 Gauge	269-**	0.095"	0.078"
14 Gauge	270-**	0.083"	0.062"
15 Gauge	271-**	0.078"	0.057"
16 Gauge	272-**	0.066"	0.050"
17 Gauge	273-**	0.058"	0.040"
18 Gauge	274-**	0.050"	0.031"
<b>11. Flexi-Needle Blunt End with Luer Lock and Tungsten Obturator**</b>			
14 Gauge	276-**	0.083"	0.062"
15 Gauge	277-**	0.078"	0.057"
15.5 Gauge	278-**	0.072"	0.054"
16 Gauge	278-**	0.066"	0.050"
17 Gauge	279-**	0.058"	0.040"
18 Gauge	280-**	0.050"	0.031"
<b>12. Flexi-Needle Blunt End with Luer Lock and Stainless Steel Obturator</b>			
13 Gauge	275-**	0.095"	0.078"
14 Gauge	756-**	0.083"	0.062"
15 Gauge	734-**	0.078"	0.057"
15.5 Gauge	1137-**	0.078"	0.054"
16 Gauge	750-**	0.066"	0.050"
17 Gauge	742-**	0.058"	0.040"
18 Gauge	764-**	0.050"	0.031"

**\*\* please specify length**

*continues on next page*

<i>Gauge</i>	<i>Part #</i>	<i>Outer Diameter</i>	<i>Inner Diameter</i>
<b>13. Flexi-Needle Blunt End, without Luer Lock, with Tungsten Obturator**</b>			
14 Gauge	361-**	0.083"	0.062"
15 Gauge	362-**	0.078"	0.057"
16 Gauge	363-**	0.066"	0.050"
17 Gauge	364-**	0.058"	0.040"
18 Gauge	365-**	0.050"	0.031"
<b>14. Flexi-Needle Blunt End, without Luer Lock, with Stainless Steel Obturator**</b>			
13 Gauge	360-**	0.095"	0.078"
14 Gauge	758-**	0.083"	0.062"
15 Gauge	736-**	0.078"	0.057"
16 Gauge	751-**	0.066"	0.050"
17 Gauge	744-**	0.058"	0.040"
18 Gauge	766-**	0.050"	0.031"
<b>15. Flexi-Needle Blunt End, without Luer Lock, with Rod (Tungsten Obturator without Luer Lock)**</b>			
14 Gauge	409-**	0.083"	0.062"
15 Gauge	410-**	0.078"	0.057"
15.5 Gauge	970-**	0.072"	0.054"
16 Gauge	411-**	0.066"	0.050"
17 Gauge	412-**	0.058"	0.040"
18 Gauge	413-**	0.050"	0.031"
<b>16. Flexi-Needle Blunt End without, Luer Lock, with Rod (Obturator Stainless Steel without Luer Lock)</b>			
13 Gauge	408-**	0.095"	0.078"
14 Gauge	760-**	0.083"	0.062"
15 Gauge	738-**	0.078"	0.057"
16 Gauge	754-**	0.066"	0.050"
17 Gauge	746-**	0.058"	0.040"
18 Gauge	768-**	0.050"	0.031"

**\*\* please specify length**

*continues on next page*

<i>Gauge</i>	<i>Part #</i>	<i>Outer Diameter</i>	<i>Inner Diameter</i>
<b>17. Tungsten Obturator only, with Luer Lock**</b>			
14 Gauge	295-**	0.055"	
15 Gauge	296-**	0.055"	
16 Gauge	297-**	0.042"	
17 Gauge	298-**	0.038"	
18 Gauge	299-**	0.028"	
<b>18. Stainless Steel Obturator only, with Luer Lock</b>			
13 Gauge	294-**	0.071"	
14 Gauge	732-**	0.055"	
15 Gauge	732-**	0.055"	
16 Gauge	748-**	0.042"	
17 Gauge	740-**	0.038"	
18 Gauge	762-**	0.028"	
<b>19. Rod only (Tungsten Obturator, without Luer Lock)**</b>			
14 Gauge	440-**	0.055"	
15 Gauge	440-**	0.055"	
16 Gauge	438-**	0.042"	
17 Gauge	436-**	0.038"	
18 Gauge	435-**	0.028"	
<b>20. Rod only (Stainless Steel Obturator, without Luer Lock)</b>			
13 Gauge	472-**	0.071"	
14 Gauge	731-**	0.055"	
15 Gauge	731-**	0.055"	
16 Gauge	747-**	0.042"	
17 Gauge	739-**	0.038"	
18 Gauge	761-**	0.028"	

**\*\* please specify length**

Sterile Flexi-Needles		
Part #	Needle Type	Quantity
1500	13 Gauge Flexi-Needle with Stainless Steel Obturator	Up to 5
1723	13 Gauge Printed Needle with Stainless Steel Obturator in a Pouch	Up to 5
1545	Stainless Steel Obturators for 13 Gauge Flexi Needles	Up to 5
1547	13 Gauge with Luer Lock without Obturator	Up to 10
1814	13 Gauge Blunt with Luer Lock without Obturator in a Pouch	Up to 10
1816	13 Gauge Blunt with Stainless Steel Obturator in a Pouch	Up to 5
1543	13 Gauge Sharp Needles	Up to 10
1549	14 Gauge Sharp Needles	Up to 10
1551	14 Gauge with Luer Lock without Obturator	Up to 10
1597	15 Gauge Set with Tungsten Obturator	20
1569	15 Gauge Set with Tungsten Obturator	15
900	15 Gauge Set with Tungsten Obturator	10
1036	15 Gauge Set with Tungsten Obturator	5
1605	15 Gauge Set with Stainless Steel Obturator	20
1577	15 Gauge Set with Stainless Steel Obturator	15
912	15 Gauge Set with Stainless Steel Obturator	10
1040	15 Gauge Set with Stainless Steel Obturator	5
1117	15 Gauge Set with Stainless Steel Obturator	1
1599	15 Gauge without Luer Lock with Obturator	20
1571	15 Gauge without Luer Lock with Obturator	15
903	15 Gauge without Luer Lock with Obturator	10
1037	15 Gauge without Luer Lock with Obturator	5
1634	15 Gauge without Luer Lock with Obturator in a Pouch	Up to 4
1525	15 Gauge Obturator Assembly with Tungsten Rod	Up to 5
1161	15 Gauge Obturator Assembly with Stainless Steel Rod	5 or 10
1607	15 Gauge without Luer Lock with Stainless Steel Obturator	20
1579	15 Gauge without Luer Lock with Stainless Steel Obturator	15
915	15 Gauge without Luer Lock with Stainless Steel Obturator	10
1041	15 Gauge without Luer Lock with Stainless Steel Obturator	5
1613	15 Gauge Printed without Luer Lock with Obturator	20

*continues on next page*

**Custom lengths available. Needles up to 25 cm long will be packed in a tray. Needles over 25 cm long will be packed in a pouch.**



Part #	Needle Type	Quantity
1585	15 Gauge Printed without Luer Lock with Obturator	15
940	15 Gauge Printed without Luer Lock with Obturator	10
942	15 Gauge Printed without Luer Lock with Obturator	5
1718	15 Gauge Printed, Sharp End with Tungsten Obturator in a Pouch	1
1144	15 Gauge Flexi-Needle with Luer Lock w/o Obturator	Up to 10
1051	15 Gauge Sharp Needles	10
1164	15 Gauge Sharp Needles	5
1166	15 Gauge Sharp Needles	1
1272	15 Gauge Blunt Needles	Up to 10
1769	15 Gauge Sharp End Flexi Needles 3-colors in a Pouch	18
1778	15 Gauge Flexi Needles Blunt End with Tungsten Obturator	10
1776	15 Gauge Flexi Needles Blunt End with Tungsten Obturator	5
1800	15 Gauge Flexi Needles without Luer Lock, Blunt End, with Stainless Steel Obturator in a tray	10
1798	15 Gauge Flexi Needles without Luer Lock, Blunt End, with Stainless Steel Obturator in a tray	5
1804	15 Gauge Flexi Needles without Luer Lock, Blunt End, with Tungsten Obturator	10
1802	15 Gauge Flexi Needles without Luer Lock, Blunt End, with Tungsten Obturator	5
1855	15 Gauge Flexi Needles Blunt End without Luer Lock with Stainless Steel Obturator in a Tray/Pouch	10
1859	15 Gauge Flexi Needle Blunt End with Luer Lock with Stainless Steel Obturator with Hub and Cap in a Tray	5
1303	15 Gauge Flexi-Needle with Stainless Steel Sleeve	5
1660	0.055" Dia. Tungsten Obturator with Sleeve and Cap	Up to 5
1658	0.055" Dia. Stainless Steel Obturator with Sleeve and Cap	Up to 5
1053	15.5 Gauge Sharp Needles	10
1615	15.5 Gauge Printed without Luer Lock with Obturator	20
1587	15.5 Gauge Printed without Luer Lock with Obturator	15

*continues on next page*

**Custom lengths available. Needles up to 25 cm long will be packed in a tray. Needles over 25 cm long will be packed in a pouch.**

<i>Part #</i>	<i>Needle Type</i>	<i>Quantity</i>
944	15.5 Gauge Printed without Luer Lock with Obturator	10
946	15.5 Gauge Printed without Luer Lock with Obturator	5
1527	15.5 Gauge with Luer Lock without Obturator	Up to 10
1339	15.5 Gauge Flexi-Needle with Stainless Steel Sleeve	Max of 10
1840	15.5 Gauge Flexi Needle without Luer Lock with Obturator in a Tray	10
1838	15.5 Gauge Flexi Needle without Luer Lock with Obturator in a Tray/Pouch	5
1836	15.5 Gauge Flexi Needle without Luer Lock with Obturator in a Pouch	Max of 4
1842	15.5 Gauge Flexi Needle Blunt End without Luer Lock with Stainless Steel Obturator in a Tray/Pouch	5
1844	15.5 Gauge Flexi Needle Blunt End without Luer Lock with Stainless Steel Obturator in a Tray	10
1601	16 Gauge Set with Tungsten Obturator	20
1573	16 Gauge Set with Tungsten Obturator	15
906	16 Gauge Set with Tungsten Obturator	10
1036	16 Gauge Set with Tungsten Obturator	5
1609	16 Gauge Set with Stainless Steel Obturator	20
1581	16 Gauge Set with Stainless Steel Obturator	15
918	16 Gauge Set with Stainless Steel Obturator	10
1042	16 Gauge Set with Stainless Steel Obturator	5
1212	16 Gauge Blunt End with Stainless Steel Obturator	1
1514	16 Gauge Sharpe End with Stainless Steel Obturator	Up to 4
1531	16 Gauge with Luer Lock without Obturator in a Pouch	Up to 10
1621	16 Gauge without Luer Lock with Obturator	20
1593	16 Gauge without Luer Lock with Obturator	15
974	16 Gauge without Luer Lock with Obturator	10
1055	16 Gauge Sharp Needles	10
1350	16 Gauge Printed Sharp End Flexi-Needles	5
1617	16 Gauge Printed without Luer Lock with Obturator	20
1589	16 Gauge Printed without Luer Lock with Obturator	15
948	16 Gauge Printed without Luer Lock with Obturator	10
950	16 Gauge Printed without Luer Lock with Obturator	5

*continues on next page*

**Custom lengths available. Needles up to 25 cm long will be packed in a tray. Needles over 25 cm long will be packed in a pouch.**



<i>Part #</i>	<i>Needle Type</i>	<i>Quantity</i>
1060	16 Gauge Blunt End with Luer Lock	1
1216	16 Gauge Sharp Needles	5
1179	16 Gauge Flexi-Needle without Luer Lock with Hub Stylet	5
1529	16 Gauge Tungsten Obturators	Up to 5
1352	16 Gauge Tungsten Stylets for Needles with Hub	Up to 5
1218	16 Gauge Stainless Steel Obturator	5
1214	16 Gauge Stainless Steel Obturator	1
1861	16 Gauge Flexi Needle Blunt End with Luer Lock with Stainless Steel Obturator in a Pouch	5
1603	17 Gauge Set with Tungsten Obturator	20
1575	17 Gauge Set with Tungsten Obturator	15
909	17 Gauge Set with Tungsten Obturator	10
1669	17 Gauge Set with Tungsten Obturator in a Pouch	Up to 4
1039	17 Gauge Set with Tungsten Obturator	5
1611	17 Gauge Set with Stainless Steel Obturator	20
1583	17 Gauge Set with Stainless Steel Obturator	15
921	17 Gauge Set with Stainless Steel Obturator	10
1671	17 Gauge Set with Stainless Steel Obturator in a Pouch	Up to 4
1867	17 Gauge Flexi Needle Blunt End with Luer Lock with Stainless Steel Obturator in a Tray/Pouch	5
1533	17 Gauge Tungsten Obturators with Rod	Up to 5
1043	17 Gauge Set with Stainless Steel with Obturator in a Tray	5
1698	17 Gauge Flexi Needle without Luer Lock with Stainless Steel Obturators	10
1556	17 Gauge Flexi Needle without Luer Lock with Stainless Steel Obturators	5
1057	17 Gauge Sharp Needles in a pouch	10
1220	17 Gauge Sharp Needles	5
1535	17 Gauge with Luer Lock without Obturator	Up to 10
1619	17 Gauge Printed without Luer Lock with Obturator	20
1591	17 Gauge Printed without Luer Lock with Obturator	15


*continues on next page*

**Custom lengths available. Needles up to 25 cm long will be packed in a tray. Needles over 25 cm long will be packed in a pouch.**

<i>Part #</i>	<i>Needle Type</i>	<i>Quantity</i>
952	17 Gauge Printed without Luer Lock with Obturator	10
954	17 Gauge Printed without Luer Lock with Obturator	5
1222	17 Gauge Stainless Steel Obturator for Needle	Up to 5
1623	18 Gauge Set with Tungsten Obturator	20
1595	18 Gauge Set with Tungsten Obturator	15
1276	18 Gauge Set with Tungsten Obturator	10
1278	18 Gauge Set with Tungsten Obturator	5
1541	18 Gauge with Luer Lock without Obturators	Up to 10
1280	18 Gauge Sharp Needles	Up to 10
1774	18 Gauge Set with Tungsten Obturators	Max of 4
1537	18 Gauge Tungsten Obturators in a Pouch	Up to 5
1539	18 Gauge Stainless Steel Obturators in a Pouch	Up to 5
1157	Flexi-Needle Insertion Kit	1 Kit

**Custom lengths available. Needles up to 25 cm long will be packed in a tray. Needles over 25 cm long will be packed in a pouch.**





**Best Medical International, Inc.** as author/distributor, believes the technical information and statements contained in this catalog to be reliable and correct at the time of publication. However, Best Medical International, Inc., does not assume any liability for statements, information or recommendations based on/or made in reliance on the information contained in this publication.

Data provided herein are for informational purposes only and are not intended for clinical use of any nature. Use of these data in such a manner is the sole responsibility of the user, and the author/distributor of the information makes no claim or inference as to its clinical suitability and/or accuracy. Best Medical International, Inc. assumes no warranty or liability for use of this information or products. Designs and specifications are subject to modification and/or change(s).

© 2013 Best Medical International, Inc.

## Terms and Conditions

**Ordering:** Orders may be placed with a purchase order number and transmitted by FAX, phone or e-mail.


**Payment:** Within 30 days.

**Freight:** FOB origin, freight prepaid and charged back through invoice.

**Returns:** Unused accessories may be returned to Best Medical International, Inc. within 15 days if found unsuitable for a particular application. No return or authorization number is necessary. A note stating the reason for return, name, and address of hospital will expedite processing of any credit.

**Note:** Best Medical International, Inc. assumes no liability for consequential damages of any kind for the material when used interchangeably with products of other manufacturers/suppliers or for any direct or indirect results and consequences of its use or misuse by the customer.

Best Medical International, Inc., reserves the right to modify the dimensions of accessories/components as well as the prices.



# Are you Searching for the Greatest Healthcare Company in the World?

## Look to TeamBest®!

World population is increasing by more than 100 million per year. To meet the needs of our aging and growing population, TeamBest® companies are committed to making quality healthcare affordable and accessible globally.

**Ready to Serve Globally for All Your Needs**

Arplay Medical  
Best Cyclotron Systems  
Best Medical Africa  
Best Medical Asia  
Best Medical Canada  
Best Medical China  
Best Medical Engineering  
Best Medical Europe  
Best Medical International  
Best Medical Latin America  
Best Medical Middle East  
Best Medical Parts  
Best Medical Research  
Best Medical Russia  
Best NDT  
Best NOMOS  
Best Particle Therapy  
Best Theratronics  
Best Vascular (Novoste)  
CNMC Company  
Genetherapy Pharmaceuticals  
Huestis Diagnostics  
Huestis Medical and more...

**For more information, visit [www.teambest.com](http://www.teambest.com)**

© 2013 Best Medical International, Inc.



Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA  
703 451 BEST (2378) 800 336 4970 [www.bestmedical.com](http://www.bestmedical.com) [www.teambest.com](http://www.teambest.com)

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA