





Best[®] medical *international*

Stainless Steel Needles



**Best Medical can
provide any standard
or customized
configuration made
to your exact
specifications, sterile
or non-sterile,
shipped within 24
hours — worldwide!**

Contact us:

**Best Medical International
7643 Fullerton Road
Springfield, Virginia 22153 USA**

Tel: 703-451-2378 • 800-336-4970

Fax: 703-451-5228

www.bestmedical.com • www.teambest.com



Mission Statement

Everyone deserves the best health care. Our goal is to work with medical professionals to provide the Best products and services. As a leader in innovative technology, our mission is to uphold our reputation for excellence in the health care field by developing, manufacturing, and delivering cost-effective, high quality products to benefit patients around the world.

Corporate Values

Fulfilling our mission requires a dedicated employee team—one that is committed to uncompromised quality and safety, and fostering a sense of pride in the company.

Best employees will:

- *Honor the dignity and uniqueness of each person, including fellow employees, customers, distributors, and vendors, through their actions, words, and behavior. They will treat others with respect, as they would wish to be treated.*
- *Embrace the spirit of cooperation and teamwork for which Best Medical is known, an environment which supports and recognizes effort and results.*
- *Promote and practice open and honest communication at all times.*
- *Listen to others and remain open to suggestions at all times, as we focus on achieving excellence in our work.*
- *Attract and retain leaders with exceptional abilities and integrity.*
- *Ensure all resources are managed in a responsible and effective manner.*

Best ADVANTAGE

- Best Medical International is the ONLY supplier of Iridium-192 and Gold-198 seeds.
- We are the first and only provider of Iridium-192 in nylon ribbons in a wide range of activities.
- We are the first and only company to provide custom service for Gold-198 permanent implants, as well as a wide range of activities.
- We are the first and only company that will call our customers the morning after shipping to inform them of mode of transportation and expected arrival time.
- We can supply Ir-192, Au-198, Pd-103, I-125 and Brachytherapy accessories on a single Purchase Order.
- We ship most materials within 24 hours of receipt of the order.

Catalog 2013 Table of Contents

Implant Needles 1

Curved Needles 2

Template Needles with Funnel and Cap..... 3

Template Needles..... 4

Needle Implanters..... 5


Diddler-Needle Adjustment Tool 6



Implant Needles


A Bevel

17G
Available in 10 cm, 15 cm, and 20 cm




B Bevel Sterile Part

17G 1076
14G 1333
Available in 10 cm, 15 cm, 20 cm, and 25 cm



C Bevel Sterile Part

17G 1458
14G
Available in 10 cm, 15 cm, 20 cm, and 25 cm
* 25 cm length only available in 17G



The Implant Needles are made of Stainless Steel and are hollow with different types of bevel at one end. Various lengths in different gauges are available, though the most commonly used needle is the 17G. Different types of Implant Needles are listed below:

	<i>A Bevel Part #</i>	<i>B Bevel Part #</i>	<i>C Bevel Part #</i>
17G 10 cm	A-203-10	B-202-10	C-201-10
17G 15 cm	A-203-15	B-202-15	C-201-15
17G 20 cm	A-203-20	B-202-20	C-201-20
17G 25 cm		B-202-25	C-201-25
14G 10 cm		B-200-10	C-199-10
14G 15 cm		B-200-15	C-199-15
14G 20 cm		B-200-20	C-199-20
Custom length is available upon request.			


CAUTION: This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.

Curved Needles

These are curved Stainless Steel Implant Needles with handle. The outer diameter (O.D.) is 0.095" and the inner diameter (I.D.) is 0.077". Two sizes are available: short and long. The short needle has a length of 3.5" and the long needle measures 7.25". Each handle is 4" long.

Also available: Thin wall curved short needles

Part # 5004-0002-Long



Part # 5004-0001-Short



CAUTION: This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.

Template Needles with Funnel and Cap



This 20 cm needle is specifically designed for use with Syed/Neblett templates or disposable templates. A funnel fixed to one end makes the ribbon loading easy and less time consuming. A special cap fits over the funnel to keep the ribbon in position.

Custom length is available upon request.

Caps are for one time use only.

CAUTION: This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.

Template Needles


Template Needles, 17G

With Collar  Closed End

Open Tip Part # 227-Length
 Sterile Open End Part # 1395
 Sterile Closed End Part # 1397

Without Collar  Tapered End

. Part # 196
 Sterile Open End Part # 1226
 Sterile Closed End Part # 1224

Without Funnel  Tapered End



Closed Tip Part # 198 Sterile Part # 1113
 Open Tip Part # 228 Sterile Part # 1115

These 17 Gauge Needles are used with the templates, which improve the homogeneity of radiation. Needles of different lengths, closed end or open end, with or without collar, with or without funnel, are available. Catalog numbers of different template needles are given below:

<i>Tapered End</i>	<i>Closed End</i>
Tapered without Collar 10 cm, 15 cm, 20 cm, 25 cm Part # 196.....Part # 1226	Closed in with Collar 10 cm, 15 cm, 20 cm, 25 cm Part # 197.....Part # 1397
Tapered with Collar 10 cm, 15 cm, 20 cm, 25 cm Part # 227.....Part # 1395	Closed in without Collar 10 cm, 15 cm, 20 cm Part # 193.....Part # 1224
Template Needle with Funnel & Tapered End 20 cm Part # 228.....Part # 1115 Stylet for the Above Part # 332.....Part # 1228	Template Needle (pencil point or sharp) with Funnel & Closed Tip 20 cm Part # 198.....Part # 1113

CAUTION: This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.

Needle Implanters

Four Inch Implanter	— 2 inch point.	Part # 1074-0016	
Six Inch Implanter	— 1 inch point.	Part # 1074-0017	
Six Inch Implanter	— 2 inch point.	Part # 1074-0019	
Six Inch Implanter	— 4 inch point.	Part # 1074-0020	
Eight Inch Implanter	— 2 inch point.	Part # 1074-0022	
Eight Inch Implanter	— 4 inch point.	Part # 1074-0023	
Ten Inch Implanter	— 1 cm point	Part # 1074-0031	
Ten Inch Implanter	— 2 inch point.	Part # 1074-0024	
Ten Inch Implanter	— 6 inch point.	Part # 1074-0025	
Ten Inch Implanter	— 8 inch point.	Part # 1074-0026	

***Other implanters also available.**

CAUTION: This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.

Diddler-Needle Adjustment Tool

Part # 5002-0010

Use the DIDDLER, Not Your Fingers!

The DIDDLER is a device to adjust the position of transperineal prostate seeding needles during the time of seed implantation. Often, fingers are needed to reach behind the transperineal template to readjust the needle position. The DIDDLER replaces the need to use fingers which, in most cases, are too big to reach a needle between the template and perineum.

The DIDDLER, similar in construction to a crochet needle, is made entirely of stainless steel and consists of a long hexagonally shaped handle with a slender tip. The tip thus allows a needle to be pushed, pulled or twisted left and right even if there is a needle right next to it.


Adjust the needle up, down and sideways without using your fingers!

Radiation Protection

The DIDDLER is especially effective when needle repositioning with "HOT" Pre-loaded Prostate Seeding Needles is required. Notably, the DIDDLER is most useful in the Mick technique when numerous needles are implanted either peripherally, in rows, or all at once.

The DIDDLER is approximately 8" (20 cm) long and made entirely of stainless steel. A high-polish finish provides for easy cleaning and flash sterilization.

The DIDDLER can be AUTOCLAVED!



Specifications:

Material	Stainless Steel
Length	8" (20 cm approximately)
Size	(6.4 mm diameter) Hexagonal Stock

CAUTION: This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.



Best Medical International, Inc. as author/distributor, believes the technical information and statements contained in this catalog to be reliable and correct at the time of publication. However, Best Medical International, Inc., does not assume any liability for statements, information or recommendations based on/or made in reliance on the information contained in this publication.

Data provided herein are for informational purposes only and are not intended for clinical use of any nature. Use of these data in such a manner is the sole responsibility of the user, and the author/distributor of the information makes no claim or inference as to its clinical suitability and/or accuracy. Best Medical International, Inc. assumes no warranty or liability for use of this information or products. Designs and specifications are subject to modification and/or change(s).

© 2012 Best Medical International, Inc.

Terms and Conditions

Ordering: Orders may be placed with a purchase order number and transmitted by FAX, phone or e-mail.


Payment: Within 30 days.

Freight: FOB origin, freight prepaid and charged back through invoice.

Returns: Unused accessories may be returned to Best Medical International, Inc. within 15 days if found unsuitable for a particular application. No return or authorization number is necessary. A note stating the reason for return, name, and address of hospital will expedite processing of any credit.

Note: Best Medical International, Inc. assumes no liability for consequential damages of any kind for the material when used interchangeably with products of other manufacturers/suppliers or for any direct or indirect results and consequences of its use or misuse by the customer.

Best Medical International, Inc., reserves the right to modify the dimensions of accessories/components as well as the prices.



Are you Searching for the Greatest Healthcare Company in the World?

Look to TeamBest®!

World population is increasing by more than 100 million per year. To meet the needs of our aging and growing population, TeamBest® companies are committed to making quality healthcare affordable and accessible globally.

Ready to Serve Globally for All Your Needs

Arplay Medical
Best Cyclotron Systems
Best Medical Africa
Best Medical Asia
Best Medical Canada
Best Medical China
Best Medical Engineering
Best Medical Europe
Best Medical International
Best Medical Latin America
Best Medical Middle East
Best Medical Parts
Best Medical Research
Best Medical Russia
Best NDT
Best NOMOS
Best Particle Therapy
Best Theratronics
Best Vascular (Novoste)
CNMC Company
Genetherapy Pharmaceuticals
Huestis Diagnostics
Huestis Medical and more...

For more information, visit www.teambest.com

© 2013 Best Medical International, Inc.



Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
703 451 BEST (2378) 800 336 4970 www.bestmedical.com www.teambest.com

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA