






# Best<sup>®</sup> medical *international*

## Templates



healthcare for everyone  
**TeamBest**<sup>®</sup>  
Your True Partner

**Best Medical can  
provide any standard  
or customized  
configuration made  
to your exact  
specifications, sterile  
or non-sterile,  
shipped within 24  
hours — worldwide!**

## **Contact us:**

**Best Medical International  
7643 Fullerton Road  
Springfield, Virginia 22153 USA**

**Tel: 703-451-2378 • 800-336-4970**

**Fax: 703-451-5228**

**[www.bestmedical.com](http://www.bestmedical.com) • [www.teambest.com](http://www.teambest.com)**



## Mission Statement

*Everyone deserves the best health care. Our goal is to work with medical professionals to provide the Best products and services. As a leader in innovative technology, our mission is to uphold our reputation for excellence in the health care field by developing, manufacturing, and delivering cost-effective, high quality products to benefit patients around the world.*

## Corporate Values

*Fulfilling our mission requires a dedicated employee team—one that is committed to uncompromised quality and safety, and fostering a sense of pride in the company.*

## Best employees will:

- *Honor the dignity and uniqueness of each person, including fellow employees, customers, distributors, and vendors, through their actions, words, and behavior. They will treat others with respect, as they would wish to be treated.*
- *Embrace the spirit of cooperation and teamwork for which Best Medical is known, an environment which supports and recognizes effort and results.*
- *Promote and practice open and honest communication at all times.*
- *Listen to others and remain open to suggestions at all times, as we focus on achieving excellence in our work.*
- *Attract and retain leaders with exceptional abilities and integrity.*
- *Ensure all resources are managed in a responsible and effective manner.*

# Best ADVANTAGE

- Best Medical International is the **ONLY** supplier of Iridium-192 and Gold-198 seeds.
- We are the first and only provider of Iridium-192 in nylon ribbons in a wide range of activities.
- We are the first and only company to provide custom service for Gold-198 permanent implants, as well as a wide range of activities.
- We are the first and only company that will call our customers the morning after shipping to inform them of mode of transportation and expected arrival time.
- We can supply Ir-192, Au-198, Pd-103, I-125 and Brachytherapy accessories on a single Purchase Order.
- We ship most materials within 24 hours of receipt of the order.

# Catalog 2013 Table of Contents

Syed/Neblett Gyn Template #1/Kit. . . . . 1

Syed/Neblett Gyn Template #2/Kit. . . . . 2

Syed/Neblett Gyn Template #3/Kit. . . . . 3

Syed/Neblett Prostate Template/Kit. . . . . 4

Syed/Neblett Rectal Template/Kit . . . . . 5

Syed/Neblett Urethral Template/Kit . . . . . 6

Disposable Gyn Template/Kit. . . . . 7



Disposable Prostate Template (with Needle Hole) . . . . . 9

Disposable Prostate Template (with Small Hole) . . . . . 10

Disposable Rectal Template . . . . . 11







**Hole Pattern**

**17G Syed/Neblett Gyn Template #1**

*Part # 1246-1 (sterile )*

*Part # 214-1 (non-sterile)*

**15G Syed/Neblett Gyn Template #1**

*Part # 1246-3 (sterile)*

*Part # 214-3 (non-sterile)*

**With:**

1	Template
1	15 cm Central Rod
2	Wrenches
1	Stainless Steel Ruler
1	Magnification Ring
	Extra 'O' Rings, Screws and Cap

**Syed/Neblett GynTemplate Kit #1**

*Part # 1246-2 (sterile)*

*Part # 214-2 (non-sterile)*

**Kit contains:**

1	Template
1	15 cm Central Rod
36	Needles with Funnels
5	Needles without Funnels
1	Magnification Ring
2	Wrenches
10	Stylets
1	Stainless Steel Ruler
	Extra 'O' Rings, Screws and Cap
	<b>Optional:</b> Storage Tray

**NOTE:** Specify hole size when ordering (15 gauge or 17 gauge).

**CAUTION:** This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.

## 17G Syed/Neblett Gyn Template #2

Part # 1248-1 (sterile)

Part # 215-1 (non-sterile)


## 15G Syed/Neblett Gyn Template #2

Part # 1248-3 (sterile)

Part # 215-3 (non-sterile)

### With:

1	Template
1	15 cm Central Rod
2	Wrenches
1	Stainless Steel Ruler
1	Magnification Ring
	Extra 'O' Rings, Screws and Cap




## 17G Syed/Neblett Gyn Template Kit #2

Part # 1248-2 (sterile)

Part # 215-2 (non-sterile)

### Kit contains:


1	Template
1	15 cm Central Rod
45	Needles with Funnels
5	Needles without Funnels
1	Magnification Ring
2	Wrenches
10	Stylets
1	Stainless Steel Ruler
	Extra 'O' Rings, Screws and Cap
	<b>Optional:</b> Storage Tray



**Hole Pattern**

**NOTE:** Specify hole size when ordering (15 gauge or 17 gauge).

**CAUTION:** This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.



**17G Syed/Neblett Gyn Template #3**

*Part # 1250-1 (sterile)*

*Part # 216-1 (non-sterile)*


**15G Syed/Neblett Gyn Template #3**

*Part # 1250-3 (sterile)*

*Part # 216-3 (non-sterile)*

**With:**

1	Template
1	15 cm Central Rod
2	Wrenches
1	Stainless Steel Ruler
1	Magnification Ring
	Extra 'O' Rings, Screws and Cap



**Hole Pattern**

**17G Syed/Neblett Gyn Template Kit #3**

*Part # 1250-2 (sterile)*

*Part # 216-2 (non-sterile)*

**Kit contains:**

1	Template
1	15 cm Central Rod
55	Needles with Funnels
5	Needles without Funnels
1	Magnification Ring
2	Wrenches
10	Stylets
1	Stainless Steel Ruler
	Extra 'O' Rings, Screws and Cap

**NOTE:** Specify hole size when ordering (15 gauge or 17 gauge).

**CAUTION:** This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.

## 17G Syed/Neblett Prostate Template

Part # 1260-1 (sterile)

Part # 159-1 (non-sterile)

## 15G Syed/Neblett Prostate Template

Part # 1260-3 (sterile)

Part # 159-3 (non-Sterile)

### With:

1	Template
1	Wrench
	Extra 'O' Rings and Screws



## 17G Syed/Neblett Prostate Template Kit

Part # 1260-2 (sterile)

Part # 159-2 (non-sterile)

### Kit contains:



1	Template
18	Needles with Funnels
3	Needles without Funnels
1	Wrench
5	Stylets
1	Stainless Steel Ruler
	Extra 'O' Rings, Screws and Cap
	<b>Optional:</b> Storage Tray



**Hole Pattern**

**NOTE:** Specify hole size when ordering (15 gauge or 17 gauge).

**CAUTION:** This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.



**Hole Pattern**

**17G Syed/Neblett Rectal Template**

*Part # 1256-1 (sterile)*

*Part # 236-1 (non-sterile)*

**15G Syed/Neblett Rectal Template**

*Part # 1256-3 (sterile)*

*Part # 236-3 (non-sterile)*

**With:**

1	Template
1	Wrench
	Extra 'O' Rings and Screws

**Syed/Neblett Rectal Template Kit**

*Part # 1256-2 (sterile)*

*Part # 236-2 (non-sterile)*

**Kit contains:**

1	Template
36	Needles with Funnels
3	Needles without Funnels
1	Wrench
1	Stainless Steel Ruler
	Extra 'O' Rings, Screws and Cap

**NOTE:** Specify hole size when ordering (15 gauge or 17 gauge).

**CAUTION:** This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.

## 17G Syed/Neblett Urethral Template

Part # 1264-1 (sterile)

Part # 237-1 (non-sterile)


## 15G Syed/Neblett Urethral Template

Part # 1264-3 (Sterile)

Part # 237-3 (Non-sterile)

### With:

1	Template
1	Wrench
	Extra 'O' Rings, Screws and Cap




## 17G Syed/Neblett Urethral Template Kit

Part # 1264-2 (sterile)

Part # 237-2 (non-sterile)

### Kit contains:

1	Template
18	Needles with Funnels
2	Needles without Funnels
1	Wrench
5	Stylets
1	Stainless Steel Ruler
	Extra 'O' Rings, Screws and Cap
	<b>Optional:</b> Storage Tray



**Hole  
Pattern**

**NOTE:** Specify hole size when ordering (15 gauge or 17 gauge).

**CAUTION:** This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.



<b>17G Disposable GYN Template Kit</b>	
<i>Part # 350-1 (non-sterile)</i>	
<b>Kit contains:</b>	
1	Template
1	15 cm Central Rod
1	Wrench
1	Sleeve
1	'O' Ring

<b>Sterile Disposable GYN Template Kit</b>	
<i>Part # 350-2 Green in a Pouch</i>	
<b>Kit contains:</b>	
290-2	Template
779	Accessories

<b>Non-sterile Template Only</b>	
Part #	
778	17G
474	15G

<b>Sterile Template Kit 15G</b>	
<i>Part # 790</i>	
<b>Kit contains:</b>	
788	Template
789	Accessories

**NOTE:** Specify hole size when ordering (15 gauge or 17 gauge).

**CAUTION:** This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.

Sterile GYN Template	
<b>Part #</b>	
290-2	Green in a Pouch

Non-Sterile Disposable GYN Purple Template Kit	
<b>Part #</b>	
691	1.9 mm

17G Sterile Accessories Template	
<b>Part #</b>	
779	In a Pouch

15G Sterile Disposable GYN Purple Template Kit	
<b>Part #</b>	<b>Kit contains:</b>
788	Template In a Pouch
789	Accessories in a Pouch
790	Complete Kit

15G Sterile Accessories Template	
<b>Part #</b>	
789	In a Pouch

16G Sterile Disposable GYN Template	
<b>Part #</b>	
1474-2	Green

16G Sterile Accessories Template	
<b>Part #</b>	
1476	In a Pouch



16G Sterile Disposable GYN Template Kit	
<b>Part #</b>	
1478-2	Green

15.5G Disposable GYN Template	
<b>Part #</b>	
1736	In a Pouch

15.5G Sterilized Disposable GYN Template Kit	
<b>Part #</b>	
1738	Complete Kit

**NOTE:** Specify hole size when ordering (15 gauge or 17 gauge).

**CAUTION:** This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.



### 17G Disposable Prostate Template Kit (with needle hole)

Part #	
223-1	Non-sterile
1184	Sterile

**NEEDLES:** 17 gauge flexi-needles or 17 gauge stainless steel needles could be used.

If you find difficulty in pushing the needles through the holes, dip needle in appropriate lubricant and use.

### 15G Disposable Prostate Template (needle hole in a pouch)

Part #	
1740	Sterile
223-4	Non-sterile

**NEEDLES:** Only 15 gauge flexi-needles or can be used with this template. If you find difficulty in pushing the needles through the holes, dip needle in appropriate lubricant and use.

This disposable prostate template is made of Medical Grade Silicone Elastomer. The hole pattern is as shown in the figure.

**This disposable template is for one time use only.**

**CAUTION:** This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.

## 17G Disposable Prostate Template (with small hole)

Part #	
473-1	Non-sterile
1186	Sterile

**NEEDLES:** 17 gauge flexi-needles or 17 gauge stainless steel needles could be used.

If you find difficulty in pushing the needles through the holes, dip needle in appropriate lubricant and use.

## 15G Prostate Template

Part #	
776-2	Non-sterile
726-2	Sterile

**NEEDLES:** 15 gauge flexi-needles or 15 gauge stainless steel needles could be used.


If you find difficulty in pushing the needles through the holes, dip needle in appropriate lubricant and use.

## 0.0555" Needle Hole Prostate Template Blue

Part #	
1853	Non-sterile
1857	Sterile

**NEEDLES:** 0.0555" flexi-needles or 0.0555" gauge stainless steel needles could be used.

If you find difficulty in pushing the needles through the holes, dip needle in appropriate lubricant and use.




This disposable prostate template is made of Medical Grade Silicone Elastomer. The hole pattern is as shown in the figure.

**This disposable template is for one time use only.**



**CAUTION:** This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.

Needle Hole




291-1

**Disposable Rectal Templates  
(with needle hole)**

Part #	
291-1	17 Gauge Non-sterile
1188	17 Gauge Sterile

Small Hole



292-1

**Disposable Rectal Templates  
(with small hole)**

Part #	
292-1	17 Gauge Non-sterile
1190	17 Gauge Sterile, Small Hole

This disposable prostate template is made of Medical Grade Silicone Elastomer.

The hole pattern is as shown in the figure.

**This disposable template is for one time use only.**



**CAUTION:** This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.

## Non-sterile Disposable Rectal Template

Part #	
345-1	16 Gauge
345-2	17 Gauge
345-3	15 Gauge
345-5	15.5 Gauge

## Sterile Disposable Rectal Template


Part #	
1084-1	17 Gauge
1084-2	15 Gauge
1084-3	16 Gauge
1084-4	15.5 Gauge

## Sterile Disposable Rectal Template Kit

Part #	
1558	15 Gauge
1678	15.5 Gauge
1664	16 Gauge
1384	17 Gauge

## 15.5G Disposable Template Accessories

Part #	
1676	Sterile in a Pouch




This disposable prostate template is made of Medical Grade Silicone Elastomer. The hole pattern is as shown in the figure.

**This disposable template is for one time use only.**

**NOTE:** Specify hole size when ordering (15 gauge or 17 gauge).

**CAUTION:** This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.




18G Templates (B&K)	
Part #	
1829	Non-sterile
1830	Sterile in a Pouch

**NOTE:** Specify hole size when ordering (15 gauge or 17 gauge).

**CAUTION:** This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.





**Best Medical International, Inc.** as author/distributor, believes the technical information and statements contained in this catalog to be reliable and correct at the time of publication. However, Best Medical International, Inc., does not assume any liability for statements, information or recommendations based on/or made in reliance on the information contained in this publication.

Data provided herein are for informational purposes only and are not intended for clinical use of any nature. Use of these data in such a manner is the sole responsibility of the user, and the author/distributor of the information makes no claim or inference as to its clinical suitability and/or accuracy. Best Medical International, Inc. assumes no warranty or liability for use of this information or products. Designs and specifications are subject to modification and/or change(s).

© 2013 Best Medical International, Inc.

## Terms and Conditions

**Ordering:** Orders may be placed with a purchase order number and transmitted by FAX, phone or e-mail.


**Payment:** Within 30 days.

**Freight:** FOB origin, freight prepaid and charged back through invoice.

**Returns:** Unused accessories may be returned to Best Medical International, Inc. within 15 days if found unsuitable for a particular application. No return or authorization number is necessary. A note stating the reason for return, name, and address of hospital will expedite processing of any credit.

**Note:** Best Medical International, Inc. assumes no liability for consequential damages of any kind for the material when used interchangeably with products of other manufacturers/suppliers or for any direct or indirect results and consequences of its use or misuse by the customer.

Best Medical International, Inc., reserves the right to modify the dimensions of accessories/components as well as the prices.



# Are you Searching for the Greatest Healthcare Company in the World?

**Look to TeamBest®!**

World population is increasing by more than 100 million per year. To meet the needs of our aging and growing population, TeamBest® companies are committed to making quality healthcare affordable and accessible globally.

**Ready to Serve Globally for All Your Needs**

Arplay Medical  
Best Cyclotron Systems  
Best Medical Africa  
Best Medical Asia  
Best Medical Canada  
Best Medical China  
Best Medical Engineering  
Best Medical Europe  
Best Medical International  
Best Medical Latin America  
Best Medical Middle East  
Best Medical Parts  
Best Medical Research  
Best Medical Russia  
Best NDT  
Best NOMOS  
Best Particle Therapy  
Best Theratronics  
Best Vascular (Novoste)  
CNMC Company  
Genetherapy Pharmaceuticals  
Huestis Diagnostics  
Huestis Medical and more...

**For more information, visit [www.teambest.com](http://www.teambest.com)**

© 2013 Best Medical International, Inc.



Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA  
703 451 BEST (2378) 800 336 4970 [www.bestmedical.com](http://www.bestmedical.com) [www.teambest.com](http://www.teambest.com)

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA