



Best[™] Integrated Brachytherapy Solutions

Best[™] Syed Templates

Best Medical offers a wide variety of disposable templates for HDR and LDR applications. Specify hole size when ordering (15, 16 or 17 gauge). These products should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.




LDR Prostate Template
(Disposable)




HDR Prostate Template
(Disposable)



Rectal Template
(Disposable)




Modified Rectal Template
(Disposable)



Grid Template for B-K type
Stepper/Stabilizer*
(Disposable)



GYN Template
(Reusable)



GYN Template
(Disposable)




Modified GYN Template
(Disposable)

Best[™] Syed Multi-Channel Vaginal Applicator

This novel HDR applicator consists of a disposable central insert and an outer sleeve along with a positioning collar. Source channels travel fully to the distal end, curving around to give more dose to the apex area, allowing for optimum dose delivery.

Additionally, the outer sleeve provides a 5 mm wall, to reduce dose to the vaginal mucosa, with a 2 mm wall at the distal end of the applicator to allow for greater dose delivery. The applicator is 18 cm in length with diameters ranging from 2.0 cm to 3.5 cm.**

*BMI provides grid templates for all makes and models of steppers/stabilizers. **All dimensions are approximate.



Sonolis[®] Transducers

We offer probes for the most common types of procedures and imaging, all with the enhanced image quality and high resolution that you would expect from the **BEST** name in ultrasounds!



- 8L2A Linear Array
- 12L5A Linear Array
- 14L3 Linear Array
- 15LW4 Linear Array
- 15LA Linear Array
- 15L4A Linear Array
- 16L5 Linear Array
- 8V3 Phased Array
- 4V2A Phased Array
- 5C2A Curved Array
- 9MC3 Curved Array
- 8EC4A Endocavity
- XY-BI-Plane Phased Array
- 10EC4A Endocavity
- 10BP4 Bi-Plane
- 8BP4 Bi-Plane
- 8TE3 Trans-esophageal
- Pedoff
- 16HL7 High Frequency Linear Array