



Best™ Integrated Brachytherapy Solutions to Improve Clinical Outcomes



The Best™ Cyber Sonalis® Ultrasound Imaging System includes our exclusive SimulView™ Technology, providing simultaneous Bi-Plane views.



Best™ HDR Remote Afterloader




Best™ Treatment Planning System



Best™ Universal Chair/Table



Best™ Localization Marker Cable



Best™ Esophageal Brachytherapy Applicator



BestMick Applicator and 10/15/20/25 Cartridges



Best™ Double-Balloon Brachytherapy Applicator for Breast and Body Cavities



Huestis Fixed Lasers—Red or Green Diode



Huestis CVO-2000 Warming Oven for Thermoplastics



mobileMOSFET Wireless Dosimetry System

Portable Dosimetry System




Vascular Brachytherapy for In-stent Restenosis

VBV Designed for the Cath Lab




- Convenient and Minimal Dosimetry Calculations
 - Portable and Reusable
 - Short Treatment Time
 - Compatible with 6F Guide and 7F Guide-Extending Catheters
- Available in 30 mm, 40 mm & 60 mm Source Train Lengths



Best[™] Integrated Radiation Therapy Solutions




Best[™] Flexi Needles



Best[™] Stepper/Stabilizer



Best[™] Brachytherapy Kit



* Half Moon Buttons are available in either radiopaque (pictured above) or clear nylon.

**Catheter tubes are available in purple, green, yellow, clear or blue with either radiopaque or clear nylon buttons.


BDS Whole Body Dosimeters



BDS Extremity Ring Dosimeters



Best[™] X-ray Marker Cables



Best[™] ABT Molecular Imaging



Best[™] B35ADP

Alpha/Deuteron/Proton Cyclotron

- **Proton Particle Beam:**
1000 μ A Beam Current
up to 35 MeV Energy
- **Deuteron Particle Beam:**
500 μ A Beam Current up
to 15 MeV Energy
- **Alpha Particle Beam:**
200 μ A Beam Current
up to 35 MeV Energy